

ULK2 Mouse mAb

Catalog No: #64387

Package Size: #64387-1 50ul #64387-2 100ul


Description

Product Name	ULK2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Purified recombinant fragment of human ULK2 (AA: 1-155) expressed in E. Coli.
Target Name	ULK2
Other Names	ATG1B; Unc51.2
Accession No.	Q8IYT8
Calculated MW	112.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using ULK2 mouse mAb against NIH/3T3 (1), HepG2 (2), SK-Hep-1 (3), SK-OV-3 (4), C6 (5), PC-12 (6), and MCF-7 (7) cell lysate.



Flow cytometric analysis of HeLa cells using ULK2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.