

USP11 Mouse mAb

Catalog No: #64391

Package Size: #64391-1 50ul #64391-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	USP11 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human USP11 expressed in E. Coli.
Target Name	USP11
Other Names	UHX1
Accession No.	P51784
Calculated MW	109.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400

Images



Western blot analysis using USP11 mouse mAb against Hela (1), LNcap (2), HepG2 (3), K562 (4), A549 (5) and Jurkat (6) cell lysate.



Flow cytometric analysis of Jurkat cells using USP11 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.