

VCAM1 Mouse mAb

Catalog No: #64397

Package Size: #64397-1 50ul #64397-2 100ul


Description

Product Name	VCAM1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human VCAM1 expressed in E. Coli.
Target Name	VCAM1
Other Names	CD106; MGC99561; INCAM-100; DKFZp779G2333; VCAM1
Accession No.	P19320
Calculated MW	81kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using VCAM1 mouse mAb against HUVEC (1) and EC (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human lung cancer (A) and colon cancer (B) using VCAM1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human Placenta tissues using VCAM1 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.