

VIL1 Mouse mAb

Catalog No: #64399

Package Size: #64399-1 50ul #64399-2 100ul


Description

Product Name	VIL1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human VIL1 (AA: 1-209) expressed in E. Coli.
Target Name	VIL1
Other Names	VIL; D2S1471
Accession No.	P09327
Calculated MW	92.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using VIL1 mouse mAb against SW620 cell lysate.



Immunofluorescence analysis of HeLa cells using VIL1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of SW620 cells using VIL1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using VIL1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.