

VISTA Mouse mAb

Catalog No: #64402

Package Size: #64402-1 50ul #64402-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VISTA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human VISTA (AA: extra 33-194) expressed in HEK293 cells.
Target Name	VISTA
Other Names	VSIR; B7H5; GI24; B7-H5; PD-1H; SISP1; PP2135; C10orf54; DD1alpha
Accession No.	Q9H7M9
Calculated MW	34kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using VISTA mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.