

VP2 Mouse mAb

Catalog No: #64403

Package Size: #64403-1 50ul #64403-2 100ul


Description

Product Name	VP2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human VP2 (AA: 296-438) expressed in E. Coli.
Target Name	VP2
Calculated MW	61kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using VP2 mouse mAb against A431 (1) and BCBL-1 (2) cell lysate.



Immunofluorescence analysis of HeLa cells using VP2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using VP2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using VP2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.