

WAS Mouse mAb

Catalog No: #64405

Package Size: #64405-1 50ul #64405-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	WAS Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human WAS (AA: 57-170) expressed in E. Coli.
Target Name	WAS
Other Names	THC; IMD2; SCN1; THC1; WASP
Accession No.	P42768
Calculated MW	53kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of Hela cells using WAS mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using WAS mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using WAS mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.