

WDR5 Mouse mAb

Catalog No: #64407

Package Size: #64407-1 50ul #64407-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	WDR5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human WDR5 expressed in E. Coli.
Target Name	WDR5
Other Names	SWD3; BIG-3; WDR5
Accession No.	P61964
Calculated MW	36.6kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500 - 1/2000 ELISA: 1/10000

Images



Western blot analysis using WDR5 mouse mAb against HeLa (1) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.