

WDR66 Mouse mAb

Catalog No: #64408

Package Size: #64408-1 50ul #64408-2 100ul


Description

Product Name	WDR66 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human WDR66 (AA: 1-250) expressed in E. Coli.
Target Name	WDR66
Other Names	WDR66
Accession No.	Q8TBY9
Calculated MW	130kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HepG2 cells using WDR66 mouse mAb (green) and negative control (purple).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.