

WIF1 Mouse mAb

Catalog No: #64410

Package Size: #64410-1 50ul #64410-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	WIF1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human WIF1 expressed in E. Coli.
Target Name	WIF1
Other Names	WIF-1; WIF
Accession No.	Q9Y5W5
Calculated MW	42kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using WIF1 mouse mAb against HeLa (1), NIH/3T3 (2) and NTERA-2 (3) cell lysate.



Immunofluorescence analysis of HeLa cells using WIF1 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Immunohistochemical analysis of paraffin-embedded human ovary tumour tissues (left) and lung cancer (right) using WIF1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.