

WNT1 Mouse mAb

Catalog No: #64411

Package Size: #64411-1 50ul #64411-2 100ul


Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | WNT1 Mouse mAb                                                                           |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | WB;IHC;IF;FC                                                                             |
| Species Reactivity    | Human,Mouse                                                                              |
| Immunogen Description | Purified recombinant fragment of WNT1 expressed in E. Coli.                              |
| Target Name           | WNT1                                                                                     |
| Other Names           | INT1                                                                                     |
| Accession No.         | P04628                                                                                   |
| Calculated MW         | 41kDa                                                                                    |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide.                                        |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000

Images



Western blot analysis using WNT1 mouse mAb against NIH/3T3 (1), 3T3L1 (2) and Hela (3) cell lysate.






Figure3: Confocal Immunofluorescence analysis of Hela (left) and 3T3-L1 (right) cells using WNT1 mouse mAb (green). Red: Actin filaments have been labeled with DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of HeLa cells using WNT1 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human Adrenal tissues using WNT1 mouse mAb

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.