

WNT10B Mouse mAb

Catalog No: #64412

Package Size: #64412-1 50ul #64412-2 100ul


Description

Product Name	WNT10B Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human WNT10B expressed in E. Coli.
Target Name	WNT10B
Other Names	SHFM6; WNT-12
Accession No.	O00744
Calculated MW	43kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using WNT10B mouse mAb against Hela cell lysate (1).



Confocal Immunofluorescence analysis of PANC-1 cells using WNT10B mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human normal stomach (A), normal liver (B), normal kidney (C) and rectum cancer tissues (D) using WNT10B mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.