

WNT5A Mouse mAb

Catalog No: #64413

Package Size: #64413-1 50ul #64413-2 100ul


Description

Product Name	WNT5A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of WNT5A expressed in E. Coli.
Target Name	WNT5A
Other Names	hWNT5A
Accession No.	P41221
Calculated MW	42.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000

Images



Confocal Immunofluorescence analysis of PC-12 cells using WNT5A mouse mAb (green), showing cytoplasmic localization. Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human lung cancer (A), thyroid cancer (B), lymph node (C) and brain (D) showing cytoplasmic and extracellular matrix localization using WNT5A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.