

YWHAQ Mouse mAb

Catalog No: #64421

Package Size: #64421-1 50ul #64421-2 100ul


Description

Product Name	YWHAQ Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human,Mouse,Rat,Monkey
Immunogen Description	Purified recombinant fragment of human YWHAQ (AA: 1-245) expressed in E. Coli.
Target Name	YWHAQ
Other Names	1C5; HS1; 14-3-3
Accession No.	P27348
Calculated MW	27.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using YWHAQ mouse mAb against Hela (1), NIH/3T3 (2), C6 (3), A549 (4), COS7 (5), PC-12 (6), and HEK293 (7) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.