

ZAP70 Mouse mAb

Catalog No: #64422

Package Size: #64422-1 50ul #64422-2 100ul


Description

Product Name	ZAP70 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ZAP70 expressed in E. Coli.
Target Name	ZAP70
Other Names	SRK; STD; TZK
Accession No.	P43403
Calculated MW	70kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using ZAP70 mouse mAb against Jurkat cell lysate (1).



Immunohistochemical analysis of paraffin-embedded human lymph tumor, showing cytoplasmic localization using ZAP70 mouse mAb with DAB staining

Note: This product is for in vitro research use only and is not intended for use in humans or animals.