

ZAP70 Mouse mAb

Catalog No: #64423

Package Size: #64423-1 50ul #64423-2 100ul


Description

Product Name	ZAP70 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ZAP70 (AA: 169-390) expressed in E. Coli.
Target Name	ZAP70
Other Names	SRK; STD; TZK; STCD; IMD48; ADMIO2; ZAP-70
Accession No.	P43403
Calculated MW	69.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using ZAP70 mouse mAb against MLOT4 (1) cell lysate.



Flow cytometric analysis of Jurkat cells using ZAP70 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.