

## ZBTB16 Mouse mAb

Catalog No: #64424

Package Size: #64424-1 50ul #64424-2 100ul


## Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | ZBTB16 Mouse mAb                                                                         |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | WB;IF                                                                                    |
| Species Reactivity    | Human                                                                                    |
| Immunogen Description | Purified recombinant fragment of human ZBTB16 expressed in E. Coli.                      |
| Target Name           | ZBTB16                                                                                   |
| Other Names           | PLZF; ZNF145; ZBTB16                                                                     |
| Accession No.         | Q05516                                                                                   |
| Calculated MW         | 74kDa                                                                                    |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide.                                        |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


## Application Details

WB:1/500 - 1/2000 ICC:1/200 - 1/1000 ELISA:1/10000

## Images



Western blot analysis using ZBTB16 mouse mAb against Hela (1) cell lysate.



Immunofluorescence analysis of Hela cells using ZBTB16 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.