

ZBTB7B Mouse mAb

Catalog No: #64425

Package Size: #64425-1 50ul #64425-2 100ul


Description

Product Name	ZBTB7B Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ZBTB7B expressed in E. Coli.
Target Name	ZBTB7B
Other Names	THPOK; ZFP67; ZBTB15; c-Krox; hcKrox; ZNF857B; DKFZp686G01254; ZBTB7B
Accession No.	O15156
Calculated MW	58kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using ZBTB7B mAb against HEK293 (1,2) cell lysate.



Immunofluorescence analysis of HeLa cells using ZBTB7B mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin



Immunohistochemical analysis of paraffin-embedded human Breast tissues using anti-ZBTB7B mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.