

ZFP42 Mouse mAb

Catalog No: #64427

Package Size: #64427-1 50ul #64427-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZFP42 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human ZFP42 (AA: 249-310) expressed in E. Coli.
Target Name	ZFP42
Other Names	REX1; ZNF754
Accession No.	Q96MM3
Calculated MW	34.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using ZFP42 mouse mAb against NIH/3T3 cell lysate.



Flow cytometric analysis of HEK293 cells using ZFP42 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using ZFP42 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using ZFP42 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.