

Recombinant Human 3-ketoacyl-CoA thiolase, peroxisomal(ACAA1)

Catalog No: #AP70003

Package Size: #AP70003-1 20ug #AP70003-2 100ug #AP70003-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human 3-ketoacyl-CoA thiolase, peroxisomal(ACAA1)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:27-331aaSequence Info:Full Length of Isoform 2
Other Names	Acetyl-CoA acyltransferase;Beta-ketothiolase;Peroxisomal 3-oxoacyl-CoA thiolase
Accession No.	P09110
Calculated MW	47.9 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	LSGAPQASAADVHVHGRRTAICRAGRGGFKDTPDELLSAVMTAVLKDVNLRPEQLGDCVGNVLQPGAGAIMARIAQFLSDIPETVPLSTVNRQCSSLQAVASIAAGGIRNGSYDIGMACGITSENVAERFGISREKQDTFALASQ QKAARAQSKGCFQAEIVPVTTHDDKGTGRSITVTQDEGIRPSTTMEGLAKLKPAFKKDGSTTAGLTVSDVDI FEINEAFASQAAYCVEKLRPLPEKVNPLGGAVALGHPLGCTGARQVITLLNELKRRGKRAYGVVSMCIGTGMG AAVFEYPGN
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

An enzyme assisted RP-RPLC approach for in-depth analysis of human liver phosphoproteome.Bian Y., Song C., Cheng K., Dong M., Wang F., Huang J., Sun D., Wang L., Ye M., Zou H.J. Proteomics 96:253-262(2014)Research Topic:Metabolism

Note: This product is for in vitro research use only and is not intended for use in humans or animals.