

Recombinant Human Acrosomal protein SP-10(ACRV1)

Catalog No: #AP70006

Package Size: #AP70006-1 20ug #AP70006-2 100ug #AP70006-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Acrosomal protein SP-10(ACRV1)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:22-265aaSequence Info:Full Length
Other Names	Acrosomal vesicle protein 1
Accession No.	P26436
Calculated MW	41.8 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	QPNELSGSIDHQTSVQQLPGEFFSLENPSDAEALYETSSGLNTLSEHGSSEHGSSKHTVAEHTSGEHAESEH ASGEPAATEHAEGEHTVGEQPSGEQPSGEHLSGEQPLSELESGEQPSDEQPSGEHGSGEQPSGEQASGEQ PSGEHASGEQASGAPISSTSTGTILNCYTCAYMNDQKGKCLRGEGETCITQNSQQCMLKKIFEGGKLQFMVQGC ENMCPSMNLFSHGTRMQIICCRNQSFCKI
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Characterization of alternatively spliced human SP-10 mRNAs.Freemerman A.J., Flickinger C.J., Herr J.C.Mol. Reprod. Dev. 41:100-108(1995)
Research Topic:Signal Transduction

Note: This product is for in vitro research use only and is not intended for use in humans or animals.