

Recombinant Human CDKN2AIP N-terminal-like protein(CDKN2AIPNL)

Catalog No: #AP70154

Package Size: #AP70154-1 20ug #AP70154-2 100ug #AP70154-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CDKN2AIP N-terminal-like protein(CDKN2AIPNL)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-116aaSequence Info:Full Length
Other Names	CDKN2A-interacting protein N-terminal-like protein
Accession No.	Q96HQ2
Calculated MW	17.2 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	MVGGEAAAIVEELVSGVRQAADFAEQFRSYSESEKQWKARMEFILRHLPDYRDPDGGSGRLDQLLSLSMWV ANHFLGCSYNKDLLDKVMEMADGIEVEDLPQFTTRSELMKKHQS
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

"N-terminal acetylome analyses and functional insights of the N-terminal acetyltransferase NatB."Van Damme P., Lasa M., Polevoda B., Gazquez C., Elosegui-Artola A., Kim D.S., De Juan-Pardo E., Demeyer K., Hole K., Larrea E., Timmerman E., Prieto J., Arnesen T., Sherman F., Gevaert K., Aldabe R.Proc. Natl. Acad. Sci. U.S.A. 109:12449-12454(2012)Research Topic:Cell Biology

Note: This product is for in vitro research use only and is not intended for use in humans or animals.