

Recombinant Chicken Collagen alpha-1(XII) chain(COL12A1) ,partial

Catalog No: #AP70180

Package Size: #AP70180-1 20ug #AP70180-2 100ug #AP70180-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Chicken Collagen alpha-1(XII) chain(COL12A1) ,partial
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:139-311aaSequence Info:Partial
Other Names	Fibrochimerin
Accession No.	P13944
Calculated MW	23.3 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	DLVFLVDGWSVGRNNFRYILDFMVALVSAFDIGEEKTRVGVVQYSSDTRTEFNLNQYFRRSDLLDAIKRIPYK GGNTMTGEAIDYLVKNTFTESAGARKGFPKVAIVITDGKAQDEVEIPARELRNIGVEVFSLGIIKAADAKELKLIAS QPSLKHVFNVANFDGIVDIQNEI
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Type XII collagen interacts with type I collagen-containing fibrils, the COL1 domain could be associated with the surface of the fibrils, and the COL2 and NC3 domains may be localized in the pericellular matrix.

References

Type XII collagen. A large multidomain molecule with partial homology to type IX collagen.Gordon M.K., Gerecke D.R., Dublet B., van der Rest M., Olsen B.R.J. Biol. Chem. 264:19772-19778(1989)Research Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.