

Recombinant Human Cancer, testis antigen 2(CTAG2)

Catalog No: #AP70203

Package Size: #AP70203-1 20ug #AP70203-2 100ug #AP70203-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cancer, testis antigen 2(CTAG2)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-210aaSequence Info:Full Length
Other Names	Autoimmunogenic cancer, testis antigen NY-ESO-2; Cancer, testis antigen 6.2 ;CT6.2L antigen family member 1 ;LAGE-1
Accession No.	O75638
Calculated MW	37.1 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MQAEGQGTTGGSTGDADGPGGPGIPDGPGGNAGGPGEAGATGGRGPRGAGAARASGPRGGAPRGPHGGA ASAQDGRCPGCGARRPDSRLQLHITMPFSSPMEAELVRRILSRDAAPLPRPGAVLKDFTVSGNLLFMSVRDQ DREGAGRMRVVGWGLGSASPEGQKARDLRTPKHKVSEQRPGTPGPPPEGAQGDGCRGVAFNVMFSAPH I
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Quantitative phosphoproteomics reveals widespread full phosphorylation site occupancy during mitosis.Olsen J.V., Vermeulen M., Santamaria A., Kumar C., Miller M.L., Jensen L.J., Gnad F., Cox J., Jensen T.S., Nigg E.A., Brunak S., Mann M.Sci. Signal. 3:RA3-RA3(2010)Research Topic:Cancer

Note: This product is for in vitro research use only and is not intended for use in humans or animals.