

Recombinant Rat Guanidinoacetate N-methyltransferase(Gamt)

Catalog No: #AP70325

Package Size: #AP70325-1 20ug #AP70325-2 100ug #AP70325-3 1mg

Description

Product Name	Recombinant Rat Guanidinoacetate N-methyltransferase(Gamt)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-236aaSequence Info:Full Length
Accession No.	P10868
Calculated MW	30.4 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	MSSSAASPLFAPGEDCGPAWRAAPAAAYDTS DTHLQILGKPVMERWETPYMHSLAAAAASRGGRVLEVGFGM AIAASRVQQAPIKEHWIIECNDGVFQRLQNWALKQPHKVPLKGLWEEEAPTLPDGHFDGILYD TYPLSEETW HTHQFNFIKTHAFRLKPGGILTYCNLT SWGELMKSKYTDITAMFEETQVPALLEAGFQRENICTEVMALVPPA DCRYYAFPQMITPLVTKH
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

crystal structures of guanidinoacetate methyltransferase ternary complexes.Komoto J., Yamada T., Takata Y., Konishi K., Ogawa H., Gomi T., Fujioka M., Takusagawa F.Biochemistry 43:14385-14394(2004)Research Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.