

# Recombinant Human Gamma-glutamylcyclotransferase(GGCT)

Catalog No: #AP70334

Package Size: #AP70334-1 20ug #AP70334-2 100ug #AP70334-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Recombinant Human Gamma-glutamylcyclotransferase(GGCT)                                                                                                                                                                                                                                                                                                                                          |
| Host Species          | E.coli                                                                                                                                                                                                                                                                                                                                                                                          |
| Purification          | Greater than 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                     |
| Immunogen Description | Expression Region:1-188aaSequence Info:Full Length                                                                                                                                                                                                                                                                                                                                              |
| Other Names           | Cytochrome c-releasing factor 21                                                                                                                                                                                                                                                                                                                                                                |
| Accession No.         | O75223                                                                                                                                                                                                                                                                                                                                                                                          |
| Calculated MW         | 37 kDa                                                                                                                                                                                                                                                                                                                                                                                          |
| Tag Info              | N-terminal 6xHis-SUMO-tagged                                                                                                                                                                                                                                                                                                                                                                    |
| Target Sequence       | MANSGCKDVTGPDEESFLYFAYGSNLLTERIHLRNPSAAFFCVARLQDFKLDGNSQGKTSQTWGGIATIFQ<br>SPGDEVWGVVWKMNKSNNLSLDEQEGVKSGMYVVEVKVATQEGKEITCRSYLMTNYESAPPSPQYKKIICM<br>GAKENGLPLEYQEKLKAIEPNDDYTGKVSSEEIEDIHKGETQTL                                                                                                                                                                                              |
| Formulation           | Tris-based buffer50% glycerol                                                                                                                                                                                                                                                                                                                                                                   |
| Storage               | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br><br>Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |

## Background

Catalyzes the formation of 5-oxoproline from gamma-glutamyl dipeptides and may play a significant role in glutathione homeostasis. Induces release of cytochrome c from mitochondria with resultant induction of apoptosis.

## References

Complete sequencing and characterization of 21,243 full-length human cDNAs.Ota T., Suzuki Y., Nishikawa T., Otsuki T., Sugiyama T., Irie R., Wakamatsu A., Hayashi K., Sato H., Nagai K., Kimura K., Makita H., Sekine M., Obayashi M., Nishi T., Shibahara T., Tanaka T., Ishii S. , Yamamoto J., Saito K., Kawai Y., Isono Y., Nakamura Y., Nagahari K., Murakami K., Yasuda T., Iwayanagi T., Wagatsuma M., Shiratori A., Sudo H., Hosoiri T., Kaku Y., Kodaira H., Kondo H., Sugawara M., Takahashi M., Kanda K., Yokoi T., Furuya T., Kikkawa E., Omura Y., Abe K., Kamihara K., Katsuta N., Sato K., Tanikawa M., Yamazaki M., Ninomiya K., Ishibashi T., Yamashita H., Murakawa K., Fujimori K., Tanai H., Kimata M., Watanabe M., Hiraoka S., Chiba Y., Ishida S., Ono Y., Takiguchi S., Watanabe S., Yosida M., Hotuta T., Kusano J., Kanehori K., Takahashi-Fujii A., Hara H., Tanase T.-O., Nomura Y., Togiya S., Komai F., Hara R., Takeuchi K., Arita M., Imose N., Musashino K., Yuuki H., Oshima A., Sasaki N., Aotsuka S., Yoshikawa Y., Matsunawa H., Ichihara T., Shiohata N., Sano S., Moriya S., Momiyama H., Satoh N., Takami S., Terashima Y., Suzuki O., Nakagawa S., Senoh A., Mizoguchi H., Goto Y., Shimizu F., Wakebe H., Hishigaki H., Watanabe T., Sugiyama A., Takemoto M., Kawakami B., Yamazaki M., Watanabe K., Kumagai A., Itakura S., Fukuzumi Y., Fujimori Y., Komiyama M., Tashiro H., Tanigami A., Fujiwara T., Ono T., Yamada K., Fujii Y., Ozaki K., Hirao M., Ohmori Y., Kawabata A., Hikiji T., Kobatake N., Inagaki H., Ikema Y., Okamoto S., Okitani R., Kawakami T., Noguchi S., Itoh T., Shigeta K., Senba T., Matsumura K., Nakajima Y., Mizuno T., Morinaga M., Sasaki M., Togashi T., Oyama M., Hata H., Watanabe M., Komatsu T., Mizushima-Sugano J., Satoh T., Shirai Y., Takahashi Y., Nakagawa K., Okumura K., Nagase T., Nomura N., Kikuchi H., Masuho Y., Yamashita R., Nakai K., Yada T., Nakamura Y., Ohara O., Isogai T., Sugano S.Nat. Genet. 36:40-45(2004)Research Topic:Metabolism

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.