

Recombinant Sheep Interferon gamma(IFNG)

Catalog No: #AP70427

Package Size: #AP70427-1 20ug #AP70427-2 100ug #AP70427-3 1mg

Description

Product Name	Recombinant Sheep Interferon gamma(IFNG)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:24-166aSequence Info:Full Length
Accession No.	P17773
Calculated MW	20.9 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	QGPFFKEIENLKEYFNASNPDAKGGPLFSEILKNWKEESDKKIIQSQIVSFYFKLFENLKDQVIQRSMDDIKQD MFQKFLNGSSEKLEDFKRLIQIPVDDLQIQRKAINELIKVMNDLSPKSNLRKRKRSQNLFRGRASM
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Produced by lymphocytes activated by specific antigens or mitogens. IFN-gamma, in addition to having antiviral activity, has important immunoregulatory functions. It is a potent activator of macrophages, it has antiproliferative effects on transformed cells and it can potentiate the antiviral and antitumor effects of the type I interferons.

References

The molecular cloning of the ovine gamma-interferon cDNA using the polymerase chain reaction.McInnes C.J., Logan M., Redmond J., Entrican G., Baird G.D.Nucleic Acids Res. 18:4012-4012(1990)Research Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.