

Recombinant Macaca mulatta Uncharacterized protein(IL17A)

Catalog No: #AP70434

Package Size: #AP70434-1 20ug #AP70434-2 100ug #AP70434-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Macaca mulatta Uncharacterized protein(IL17A)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:24-155aaSequence Info:Full Length
Accession No.	F6T3G5
Calculated MW	42.1 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	GIAIPRNP GCPNSEDKTFPRTVMVNLNIHNRNTNTNPKRSSDYNNRSTSPWNLHRNEDPERYPSVIWEAKCR HLGCVNADGNVDYHMNSVPIQQEILVLRREPRHCPNSFRLEKILVSVGCTCVTPIVHHVA
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Genome sequencing and comparison of two nonhuman primate animal models, the cynomolgus and Chinese rhesus macaques.Yan G., Zhang G., Fang X., Zhang Y., Li C., Ling F., Cooper D.N., Li Q., Li Y., van Gool A.J., Du H., Chen J., Chen R., Zhang P., Huang Z., Thompson J.R., Meng Y., Bai Y. , Wang J., Zhuo M., Wang T., Huang Y., Wei L., Li J., Wang Z., Hu H., Yang P., Le L., Stenson P.D., Li B., Liu X., Ball E.V., An N., Huang Q., Zhang Y., Fan W., Zhang X., Li Y., Wang W., Katze M.G., Su B., Nielsen R., Yang H., Wang J., Wang X., Wang J.Nat. Biotechnol. 29:1019-1023(2011) Research Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.