

Recombinant Human Serine,threonine-protein phosphatase 2B catalytic subunit gamma isoform(PPP3CC),partial

Catalog No: #AP70657

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP70657-1 20ug #AP70657-2 100ug #AP70657-3 1mg

Description

Product Name	Recombinant Human Serine,threonine-protein phosphatase 2B catalytic subunit gamma isoform(PPP3CC),partial
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:2-512aaSequence Info:Partial
Other Names	CAM-PRP catalytic subunitCalcineurin, testis-specific catalytic subunitCalmodulin-dependent calcineurin A subunit gamma isoform
Accession No.	P48454
Calculated MW	62 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	SGRRFHLSTTDRVIKAVFPPTQRLTFKEVFENGKPKVDVLKNHLVKEGRLEEEVALKIINDGAAILRQEKTMI VDAPITVCGDIHGQFFDLMKLFEVGGSPSNTRYLFLGDYVDRGYFSIECVLYLWSLKINHPKTLFLLRGNHECR HLTDYFTFKQCEKIKYSEQVYDACMETFDCLPLAALLNQQLFCVHGGMSPEITSLDDIRKLDRFTEPPAFGPVC DLLWSDPSEDYGNEKTLEHYTHNTVRGCSYFYSPAVCEFLQNNNLLSIIRAHEAQDAGYRMYRKSQATGFP SLITIFSAPNYLDVYNKAAVLKYENNVNIRQFNCSPHPYWLPNFMDFVFTWSLPFVGEKVTEMLVNVNLCSD DELISDDEAEGSTTVRKEIIRNKIRAIGKMARVFSILRQESVLTGLTPTGTPLGLVLSGGKQTIETATVEAVE AREAIRGFSLQHKIRSFEEARGLDRINERMPPRKDSIHAGGPMKSVTSAHSHAAHRSDQGKKAHS
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Calcium-dependent, calmodulin-stimulated protein phosphatase. This subunit may have a role in the calmodulin activation of calcineurin.

References

DNA sequence and analysis of human chromosome 8.Nusbaum C., Mikkelsen T.S., Zody M.C., Asakawa S., Taudien S., Garber M., Kodira C.D., Schueler M.G., Shimizu A., Whittaker C.A., Chang J.L., Cuomo C.A., Dewar K., FitzGerald M.G., Yang X., Allen N.R., Anderson S., Asakawa T. , Blechschmidt K., Bloom T., Borowsky M.L., Butler J., Cook A., Corum B., DeArellano K., DeCaprio D., Dooley K.T., Dorris L. III, Engels R., Gloeckner G., Hafez N., Hagopian D.S., Hall J.L., Ishikawa S.K., Jaffe D.B., Kamat A., Kudoh J., Lehmann R., Lokitsang T., Macdonald P., Major J.E., Matthews C.D., Mauceli E., Menzel U., Mihalev A.H., Minoshima S., Murayama Y., Naylor J.W., Nicol R., Nguyen C., O'Leary S.B., O'Neill K., Parker S.C.J., Polley A., Raymond C.K., Reichwald K., Rodriguez J., Sasaki T., Schilhabel M., Siddiqui R., Smith C.L., Sneddon T.P., Talamas J.A., Tenzin P.,

Note: This product is for in vitro research use only and is not intended for use in humans or animals.