

Recombinant Chicken Syntaxin-6(STX6),partial

Catalog No: #AP70814

Package Size: #AP70814-1 20ug #AP70814-2 100ug #AP70814-3 1mg

Description

Product Name	Recombinant Chicken Syntaxin-6(STX6),partial
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:4-231aaSequence Info:Partial
Accession No.	Q5ZL19
Calculated MW	30.1 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	EDPFFVVKGEVQKAVNTAQGLFQRWTELLQDPSIATREEIDWTTNELRNNLRSIEWDLEDLDETISIVEANPRK FNLDATELGIRKSFITSTRQVVREMKDQMSNSSMQALAERKNRQALLGESSSQSWSSGPDKYSRDLRELQLA NSHFIEEQQAQQQLIVEQQDEQLELVSGSIGVLKNMSQRIGGELEEQAVMLDDFSHELDSTHSRLDNVMKKLA KVSHMTSDR
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in intracellular vesicle trafficking.

References

Full-length cDNAs from chicken bursal lymphocytes to facilitate gene function analysis.Caldwell R.B., Kierzek A.M., Arakawa H., Bezzubov Y., Zaim J., Fiedler P., Kutter S., Blagodatski A., Kostovska D., Koter M., Plachy J., Carninci P., Hayashizaki Y., Buerstedde J.-M.Genome Biol. 6:R6.1-R6.9(2005)Research Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.