

Recombinant Human Transmembrane protein 98(TM98)

Catalog No: #AP70862

Package Size: #AP70862-1 20ug #AP70862-2 100ug #AP70862-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Transmembrane protein 98(TM98)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:25-226aaSequence Info:Full Length
Other Names	Protein TADA1
Accession No.	Q9Y2Y6
Calculated MW	26.2 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	LVCRQRYCRPRDLLQRYDSKPIVDLIGAMETQSEPSELEDDVVITNPHIEAILENEDWIEDASGLMSHCIAILKI CHTLTEKLVAMTMGSGAKMKTSASVSDIIVAKRISPRVDDVVKSMYPPLDPKLLDARTTALLSVSHLVLVTR NACHLTGGLDWIDQSLSAEEHLEVLREAALASEPDKGLPGPEGFLQEQSAI
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Teramoto T.Towards a catalog of human genes and proteins sequencing and analysis of 500 novel complete protein coding human cDNAs.Wiemann S., Weil B., Wellenreuther R., Gassenhuber J., Glassl S., Ansorge W., Boecher M., Bloecker H., Bauersachs S., Blum H., Lauber J., Duesterhoeft A., Beyer A., Koehler K., Strack N., Mewes H.-W., Ottenwaelder B., Obermaier B. , Tampe J., Heubner D., Wambutt R., Korn B., Klein M., Poustka A.Genome Res. 11:422-435(2001)Research Topic:Cell Biology

Note: This product is for in vitro research use only and is not intended for use in humans or animals.