

Recombinant Human Tumor protein D52(TPD52)

Catalog No: #AP70873

Package Size: #AP70873-1 20ug #AP70873-2 100ug #AP70873-3 1mg

Description

Product Name	Recombinant Human Tumor protein D52(TPD52)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-184aaSequence Info:Full Length of Isoform 2
Other Names	Protein N8
Accession No.	P55327
Calculated MW	35.9 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MDRGEQGLLRTPVPEEGEDVAATISATETLSEEEQEELRRELAKVEEEIQTL SQVLAAKEKHLAEIKRKLGINSLQELKQNIAGWQDVTATSAYKKTSETLSQAGQKASAAFSSVGSVITKKLEDVKNSPTFKSFEEKVENLKSIVGGTKPAGGDFGEVLNSAANASATTTEPLPEKTQESL
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Transcription variants of the prostate-specific PrLZ gene and their interaction with 14-3-3 proteins.Wang R., He H., Sun X., Xu J., Marshall F.F., Zhau H., Chung L.W., Fu H., He D.Biochem. Biophys. Res. Commun. 389:455-460(2009)Research Topic:Cancer

Note: This product is for in vitro research use only and is not intended for use in humans or animals.