

Recombinant Human Protein yippee-like 3(YPEL3)

Catalog No: #AP70926

Package Size: #AP70926-1 20ug #AP70926-2 100ug #AP70926-3 1mg

Description

Product Name	Recombinant Human Protein yippee-like 3(YPEL3)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-119aaSequence Info:Full Length
Accession No.	P61236
Calculated MW	29.6 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MVRISKPKTFQAYLDDCHRRYSCAHCRAHLANHDDLISKSFQGSQGRAYLFNSVVNVGCGPAEERVLLTGLH AVADIHCENCKTTLGWKYEQAFESSQKYKEGKYIIELNHHMIKDNGWD
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in proliferation and apoptosis in myeloid precursor cells.

References

Identification and characterization of a novel gene family YPEL in a wide spectrum of eukaryotic species.Hosono K., Sasaki T., Minoshima S., Shimizu N.Gene 340:31-43(2004)Research Topic:Epigenetics and Nuclear Signaling

Note: This product is for in vitro research use only and is not intended for use in humans or animals.