

Recombinant Human Zinc finger protein 114(ZNF114),partial

Catalog No: #AP70927

Package Size: #AP70927-1 20ug #AP70927-2 100ug #AP70927-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Zinc finger protein 114(ZNF114),partial
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-414aaSequence Info:Partial
Accession No.	Q8NC26
Calculated MW	63.4 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MSQDSVTFADVAVNFTKEEWTLDDPAQRNLYRDVMLENSRNLAFIDWATPCKTKDATPQPDILPKRTFPEANR VCLTSISSQHSTLREDWRCPKTEEPHRQGVNNVKPPAFAPEKDESPVSICEDHEMRNHSKPTCRLVPSQGDS IRQCILTRDSSIFKYNPVLNDSQKTHENNEDDGLVGNWVPCGRKTELKSSWTGSGQNTVHHIRDEIDTGA NRHQRNPFKGAFREDGSLRAHNTHGREKMYDFTQCENSRNNSIHAMQMQLYTAETNKKDCQTGATSANA PNSGSHKSHCTGEKTHKCPECGRAFFYQSFLMRHMKIHTGEKPYECGKCGKAFRYSLHLNKLHRLKHVVQKK PYECEECGKVIRESSKYTHIRSHTGEKPYKCKTCGKDFAKSSGLKKHLKTHKDEK
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May be involved in transcriptional regulation.

References

Complete sequencing and characterization of 21,243 full-length human cDNAs.Ota T., Suzuki Y., Nishikawa T., Otsuki T., Sugiyama T., Irie R., Wakamatsu A., Hayashi K., Sato H., Nagai K., Kimura K., Makita H., Sekine M., Obayashi M., Nishi T., Shibahara T., Tanaka T., Ishii S. , Yamamoto J., Saito K., Kawai Y., Isono Y., Nakamura Y., Nagahari K., Murakami K., Yasuda T., Iwayanagi T., Wagatsuma M., Shiratori A., Sudo H., Hosoiri T., Kaku Y., Kodaira H., Kondo H., Sugawara M., Takahashi M., Kanda K., Yokoi T., Furuya T., Kikkawa E., Omura Y., Abe K., Kamihara K., Katsuta N., Sato K., Tanikawa M., Yamazaki M., Ninomiya K., Ishibashi T., Yamashita H., Murakawa K., Fujimori K., Tanai H., Kimata M., Watanabe M., Hiraoka S., Chiba Y., Ishida S., Ono Y., Takiguchi S., Watanabe S., Yosida M., Hotuta T., Kusano J., Kanehori K., Takahashi-Fujii A., Hara H., Tanase T.-O., Nomura Y., Togiya S., Komai F., Hara R., Takeuchi K., Arita M., Imose N., Musashino K., Yuuki H., Oshima A., Sasaki N., Aotsuka S., Yoshikawa Y., Matsunawa H., Ichihara T., Shiohata N., Sano S., Moriya S., Momiyama H., Satoh N., Takami S., Terashima Y., Suzuki O., Nakagawa S., Senoh A., Mizoguchi H., Goto Y., Shimizu F., Wakebe H., Hishigaki H., Watanabe T., Sugiyama A., Takemoto M., Kawakami B., Yamazaki M., Watanabe K., Kumagai A., Itakura S., Fukuzumi Y., Fujimori Y., Komiyama M., Tashiro H., Tanigami A., Fujiwara T., Ono T., Yamada K., Fujii Y., Ozaki K., Hirao M., Ohmori Y., Kawabata A., Hikiji T., Kobatake N., Inagaki H., Ikema Y., Okamoto S., Okitani R., Kawakami T., Noguchi S., Itoh T., Shigeta K., Senba T., Matsumura K., Nakajima Y., Mizuno T., Morinaga M., Sasaki M., Togashi T., Oyama M., Hata H., Watanabe M., Komatsu T., Mizushima-Sugano J.,

Note: This product is for in vitro research use only and is not intended for use in humans or animals.