

Recombinant Human Zinc finger protein 155(ZNF155)

Catalog No: #AP70928

Package Size: #AP70928-1 20ug #AP70928-2 100ug #AP70928-3 1mg

Description

Product Name	Recombinant Human Zinc finger protein 155(ZNF155)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-538aaSequence Info:Full Length
Accession No.	Q12901
Calculated MW	78.1 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MTTFKEAVTFKDVAVVFTEELGLLDPAQRKLYRDVMLENFRNLLSVGHQPFHQDTCHFLREEKFWMMGTAT QREGNSGGKIQTELESVPEAGAHEEWSCQIQIWEQIAKDLTRSQDSIINNSQFFENGDVPSQVEAGLPTIHTGQ KPSQGGKCKQSISDVPIFDLPQQLYSEEKSYTCDECGKSICYISALHVHQRVHVGEKLFMCDVCGKEFSQSSH LQTHQRVHTGEKPFKCEQCGKGFSSRSALNVHRKLHTGEKPYICEACGKAFIHDSQLKEHKRIHTGEKPFKCD ICGKTFYFRSRLKSHSMVHTGEKPFRCDCDKSFHQRSALNRHCMVHTGEKPYRCEQCGKGFGRDLDFYKH QVVHTGEKPYNCKEKGKSFWRSSCLLNHQRVHSGEKSFKCEEKGKGFYTNSQLSSHQRSHSGEKPYKCEE CGKGYVTNFNLDLHQRVHTGERPYNCKEKGKNSRASSILNHKRLHCQKPKFCEDCGKRLVHRTYRKDQP RDYSGENPSKCEDCGRRYKRRLNLDILLSLFLNDT
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May be involved in transcriptional regulation.

References

The consensus coding sequences of human breast and colorectal cancers.Sjoebloom T., Jones S., Wood L.D., Parsons D.W., Lin J., Barber T.D., Mandelker D., Leary R.J., Ptak J., Silliman N., Szabo S., Buckhaults P., Farrell C., Meeh P., Markowitz S.D., Willis J., Dawson D., Willson J.K.V. , Gazdar A.F., Hartigan J., Wu L., Liu C., Parmigiani G., Park B.H., Bachman K.E., Papadopoulos N., Vogelstein B., Kinzler K.W., Velculescu V.E.Science 314:268-274(2006)Research Topic:Epigenetics and Nuclear Signaling

Note: This product is for in vitro research use only and is not intended for use in humans or animals.