

Recombinant Epstein-Barr virus (strain GD1) (HHV-4) DUTPase(BLLF3)

Catalog No: #AP70988

Package Size: #AP70988-1 20ug #AP70988-2 100ug #AP70988-3 1mg

Description

Product Name	Recombinant Epstein-Barr virus (strain GD1) (HHV-4) DUTPase(BLLF3)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-278aaSequence Info:Full Length
Accession No.	K9US42
Calculated MW	46.9 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MEACPHIRYAFQNDKLLLQQASVGRLLVNKTTILLRPMKTTTVDLGLYARPPEGHGLMLWGSTSRPVTSHVG IIDPGYTGELRLILQNQRRYNSTLRPSELKIHAAAFRYATPQMEEDKGPINHPQYPGDVGLDVSLPKDLALFPHQ TVSVTLTVPPPSIPHHRPTIFGRSGLAMQGILVKPCRWRRGVDVSLTNFSDQTVFLNKYRRFCQLVYLHKHH LTSFYSPHSDAGVLGPRSLFRWASCAFEVPSLAMGDSGLSEALEGRQGRGFGSSGQ
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Identification and characterization of EBV genomes in spontaneously immortalized human peripheral blood B lymphocytes by NGS technology.Lei H., Li T., Hung G.C., Li B., Tsai S., Lo S.C.BMC Genomics 14:804-804(2013) Research Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.