

Recombinant?Schistosoma?mansoni?Major?egg?anti gen

Catalog No: #AP70991

Package Size: #AP70991-1 20ug #AP70991-2 100ug #AP70991-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant?Schistosoma?mansoni?Major?egg?antigen
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-354aaSequence Info:Full Length
Other Names	p40
Accession No.	P12812
Calculated MW	43.4 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	MSGGKQHNAVSIPVNREQRSFEKQRRDLLTGLEHGGGAHRGNSIAPYTEDWPSTVDNWIDSSWKRWDDDM RRLRRGMFALLPLDTFSIGILENPFALMHQMDRQIQDIRERMGSLDVPSTGSVNDFLKDAYEVEGEDGKVHFKV RFDAQGFAPQDINVTSSNRVTVHAKKETTTDGRKCSREFCRMVQLPKSIDDSQLKCRMTDDGVLMLLEAPVK VDQNQSLTLNESGQVAVRPKSDNQIKAVPASQALVAKGVHGLSYVDDGSGGKRLHVEVAVDVPVYKPEDLFVN VDSNRVVVSGRHHKQKSDQHGRSSSFAEFSQSYAIPETVDPLSVSAQVVGNTLVLEAPLEKQHAITH
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

"Sequence and expression of a major egg antigen from Schistosoma mansoni. Homologies to heat shock proteins and alpha-crystallins."Nene V., Dunne D.W., Johnson K.S., Taylor D.W., Cordingley J.S.Mol. Biochem. Parasitol. 21:179-188(1986) Research Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.