

Recombinant Human MAP3K7 C-terminal-like protein(MAP3K7CL)

Catalog No: #AP71096

Package Size: #AP71096-1 20ug #AP71096-2 100ug #AP71096-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MAP3K7 C-terminal-like protein(MAP3K7CL)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-142aaSequence Info:Full Length of Isoform D
Other Names	TAK1-like protein
Accession No.	P57077
Calculated MW	32.4 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MISTARVPADKPVRIAFSLNDASDDTPPEDSIPLVFPELDQQLQPLPPCHDSEESMEVFKQHCQIAEEYHEVKK EITLLEQRKKELIAKLDQAEKEKVDAAELVREFEALTEENRTLRLAQSQCVEQLEKLRIQYQKRQGSS
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

The DNA sequence of human chromosome 21.Hattori M., Fujiyama A., Taylor T.D., Watanabe H., Yada T., Park H.-S., Toyoda A., Ishii K., Totoki Y., Choi D.-K., Groner Y., Soeda E., Ohki M., Takagi T., Sakaki Y., Taudien S., Blechschmidt K., Polley A. , Menzel U., Delabar J., Kumpf K., Lehmann R., Patterson D., Reichwald K., Rump A., Schillhabel M., Schudy A., Zimmermann W., Rosenthal A., Kudoh J., Shibuya K., Kawasaki K., Asakawa S., Shintani A., Sasaki T., Nagamine K., Mitsuyama S., Antonarakis S.E., Minoshima S., Shimizu N., Nordsiek G., Hornischer K., Brandt P., Scharfe M., Schoen O., Desario A., Reichelt J., Kauer G., Bloecker H., Ramser J., Beck A., Klages S., Hennig S., Riesselmann L., Dagand E., Wehrmeyer S., Borzym K., Gardiner K., Nizetic D., Francis F., Lehrach H., Reinhardt R., Yaspo M.-L.Nature 405:311-319(2000)Research Topic:Signal Transduction

Note: This product is for in vitro research use only and is not intended for use in humans or animals.