

Recombinant Staphylococcus aureus Enterotoxin type B(entB)

Catalog No: #AP71147

Package Size: #AP71147-1 20ug #AP71147-2 100ug #AP71147-3 1mg



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Staphylococcus aureus Enterotoxin type B(entB)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:28-266aaSequence Info:Full Length of Mature Protein
Other Names	SEB
Accession No.	P01552
Calculated MW	55.4 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	ESQPDPKPDELHKSSKFTGLMENMKVLYDDNHVSAINVKSIDQFLYFDLIYSIKDTKLGNYDNVRVEFKNKDLA DKYKDKYVDVFGANYYYQCYFSKKTNDINSHQTDKRKTCMYGGVTEHNGNQLDKYRSITVRVFEDGKNLLSF DVQTNKKKVTAQELDYLTRHYLVKNKKLYEFNNSPYETGYIKFIENENSWYDMMPAPGDKFDQSKYLMMYN DNKMVDSKDVKIEVYLTTKKK
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Staphylococcal enterotoxins cause the intoxication staphylococcal food poisoning syndrome. The illness characterized by high fever, hypotension, diarrhea, shock, and in some cases death.

References

Nucleotide sequence of the enterotoxin B gene from Staphylococcus aureus.Jones C.L., Khan S.A.J. Bacteriol. 166:29-33(1986)Research
Topic:Others

Note: This product is for in vitro research use only and is not intended for use in humans or animals.