

Calnexin

Catalog No: #AP81229

Package Size: #AP81229-1 50ug #AP81229-2 100ug #AP81229-3 1mg



Description

Product Name	Calnexin
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	pp90
Accession No.	Uniprot ID:P24643
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CANX

Published Papers

Yihong Liu;Yihong Liu;Yihong Liu;Qingshan Ma;Qingshan Ma;Qingshan Ma;Muhammad Zahoor Khan;Muhammad Zahoor Khan;Muhammad Zahoor Khan;Menghan Wang;Menghan Wang;Menghan Wang;Fokun Xiang;Fokun Xiang;Fokun Xiang;Xinyue Zhang;Xinyue Zhang;Xinyue Zhang;Xiyan Kou;Xiyan Kou;Xiyan Kou;Shuhuan Li;Shuhuan Li;Shuhuan Li;Changfa Wang;Changfa Wang;Changfa Wang;Yan Li;Yan Li;Yan Li et al., Proteomic Profiling of Donkey Milk Exosomes Highlights Bioactive Proteins with Immune-Related Functions, , (2025)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.