

## Recombinant Mouse EGF Protein, C-Fc

Catalog No: #AP96077

Package Size: #AP96077-1 50ug #AP96077-2 100ug #AP96077-3 1mg

## Description

|                |                                                                                                                                                                           |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Mouse EGF Protein, C-Fc                                                                                                                                       |
| Host Species   | Mammalian Cells                                                                                                                                                           |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                           |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                           |
| Other Names    | Pro-epidermal growth factor, EGF, Epidermal growth factor, Egf                                                                                                            |
| Accession No.  | P01132                                                                                                                                                                    |
| Target Species | Mus musculus (Mouse)                                                                                                                                                      |
| Calculated MW  | 37.12 kDa                                                                                                                                                                 |
| Target Length  | Asn977-Arg1029                                                                                                                                                            |
| Tag Info       | Fc-tag                                                                                                                                                                    |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                           |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.