

Recombinant Human SH2D1A Protein, N-GST & C-His

Catalog No: #AP97487

Package Size: #AP97487-1 50ug #AP97487-2 100ug #AP97487-3 1mg

Description

Product Name	Recombinant Human SH2D1A Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	SH2 domain-containing protein 1A, Duncan disease SH2-protein, Signaling lymphocytic activation molecule-associated protein, SLAM-associated protein, T-cell signal transduction molecule SAP, SH2D1A, DSHP, SAP
Accession No.	O60880
Target Species	Homo sapiens (Human)
Calculated MW	38.88 kDa
Target Length	Asp2-Pro97
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.