

Recombinant Human GCM2 Protein, N-His

Catalog No: #AP97511

Package Size: #AP97511-1 50ug #AP97511-2 100ug #AP97511-3 1mg

Description

Product Name	Recombinant Human GCM2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Chorion-specific transcription factor GCMb, hGCMb, GCM motif protein 2, Glial cells missing homolog 2, GCM2, GCMB
Accession No.	O75603
Target Species	Homo sapiens (Human)
Calculated MW	22.74 kDa
Target Length	Met1-Tyr183
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.