

Recombinant Human CTAG2 Protein, N-GST & C-His

Catalog No: #AP97514

Package Size: #AP97514-1 50ug #AP97514-2 100ug #AP97514-3 1mg

Description

Product Name	Recombinant Human CTAG2 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Cancer/testis antigen 2, CT2, Autoimmunogenic cancer/testis antigen NY-ESO-2, Cancer/testis antigen 6.2, CT6.2, L antigen family member 1, LAGE-1, CTAG2, ESO2, LAGE1
Accession No.	O75638
Target Species	Homo sapiens (Human)
Calculated MW	34.33 kDa
Target Length	Leu87-Asp141
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.