

Recombinant Human EIF3J Protein, N-His

Catalog No: #AP97522

Package Size: #AP97522-1 50ug #AP97522-2 100ug #AP97522-3 1mg

Description

Product Name	Recombinant Human EIF3J Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Eukaryotic translation initiation factor 3 subunit J, eIF3j, Eukaryotic translation initiation factor 3 subunit 1, eIF-3-alpha, eIF3 p35, EIF3J, EIF3S1, PRO0391
Accession No.	O75822
Target Species	Homo sapiens (Human)
Calculated MW	15.50 kDa
Target Length	Leu108-Lys222
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.