

Recombinant Human TUSC2 Protein, N-GST & C-His

Catalog No: #AP97526

Package Size: #AP97526-1 50ug #AP97526-2 100ug #AP97526-3 1mg

Description

Product Name	Recombinant Human TUSC2 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Tumor suppressor candidate 2, Fusion 1 protein, Fus-1 protein, PDGFA-associated protein 2, TUSC2, C3orf11, FUS1, LGCC, PDAP2
Accession No.	O75896
Target Species	Homo sapiens (Human)
Calculated MW	34.19 kDa
Target Length	Thr46-His96
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.