

Recombinant Human PIAS2 Protein, N-His

Catalog No: #AP97528

Package Size: #AP97528-1 50ug #AP97528-2 100ug #AP97528-3 1mg

Description

Product Name	Recombinant Human PIAS2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	E3 SUMO-protein ligase PIAS2, 2.3.2.-, Androgen receptor-interacting protein 3, ARIP3, DAB2-interacting protein, DIP, E3 SUMO-protein transferase PIAS2, Msx-interacting zinc finger protein, Miz1, PIAS-NY protein, Protein inhibitor of activated STAT x, Protein inhibitor of activated STAT2, PIAS2, PIASX
Accession No.	O75928
Target Species	Homo sapiens (Human)
Calculated MW	33.39 kDa
Target Length	Leu151-Cys425
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.