

Recombinant Human TOX Protein, N-GST & C-His

Catalog No: #AP97539

Package Size: #AP97539-1 50ug #AP97539-2 100ug #AP97539-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TOX Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Thymocyte selection-associated high mobility group box protein TOX, Thymus high mobility group box protein TOX, TOX, KIAA0808
Accession No.	O94900
Target Species	Homo sapiens (Human)
Calculated MW	36.75 kDa
Target Length	Val265-Val341
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.