

Recombinant Human GLCE Protein, N-His

Catalog No: #AP97541

Package Size: #AP97541-1 50ug #AP97541-2 100ug #AP97541-3 1mg

Description

Product Name	Recombinant Human GLCE Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	D-glucuronyl C5-epimerase, 5.1.3.17, Heparan sulfate C5-epimerase, Hsepi, Heparin/heparan sulfate:glucuronic acid C5-epimerase, Heparosan-N-sulfate-glucuronate 5-epimerase, GLCE, KIAA0836
Accession No.	O94923
Target Species	Homo sapiens (Human)
Calculated MW	52.34 kDa
Target Length	Ser176-Asn617
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.