

Recombinant Human KDM4B Protein, N-His

Catalog No: #AP97543

Package Size: #AP97543-1 50ug #AP97543-2 100ug #AP97543-3 1mg

Description

Product Name	Recombinant Human KDM4B Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Lysine-specific demethylase 4B, 1.14.11.66, JmjC domain-containing histone demethylation protein 3B, Jumonji domain-containing protein 2B, [histone H3]-trimethyl-L-lysine(9) demethylase 4B, KDM4B, JHDM3B, JMJD2B, KIAA0876
Accession No.	O94953
Target Species	Homo sapiens (Human)
Calculated MW	39.54 kDa
Target Length	Met16-Asp338
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.